

AUTO CLAVE



DK-AC 010



DK-AC 010 Inner Shape



DK-AC 001 Lid


Application

■ Sterilize glassware, laboratory equipments, media and waste.

General Features

- Full Automatic System : After setting temp. and time, all procedures are done automatically.
- Ending Alarm: If all procedures are finished, there is always ending alarm.
- Special designed Silicon seal packing for lid.
- Protecting device for over temperature and over pressure.
- Adjustable valve which is controlled by hand, if necessary.

Automatic Sterilization Process



Specification

Model	DK-AC 001	DK-AC 010	DK-AC 020	DK-AC 030
Chamber Size	<i></i>	<i>⊕</i> 350 x 630 _{mm}	<i></i>	<i></i> 450 x 630 _{mm}
Capacity	45 Liter	60 Liter	80 Liter	100 Liter
Type	Full Automatic System			
Chamber Type	Round-Vertical Type			
Temp. Controller	Digital P.I.D Controller (Auto Tuning)			
Temp. Range	Ambient to 132°C			
Operating Pressure	1.2Kg/♂-121 ℃			
Timer	99min, Digital Timer			
Door	One Touch Type (Single Screw Door)			
Pressure Gauge	Mechanical Type 0-3Kg/c㎡			
Air Exhaust	Adjustable Handle Valve			
Safety Device	Over Heat Protector, Over Pressure Protector, Alarm			
Basket	Standard Type 2EA			
Power	2.5KW	3KW	4KW	5KW